



**BUILDING THE MOST RUGGED PRODUCTS WITH UNRIVALED CRAFTSMANSHIP AND CUSTOMIZATION**

## THE CROWN JEWEL OF THE BRAUN LINEUP

Trusted by departments nationwide, the Chief XL delivers the rugged reliability and precision craftsmanship that defines the Braun standard of excellence.

Available in both Type I and Type III configurations with seven chassis options in 2WD or 4WD, the Chief XL adapts to your department's specific operational requirements.

Built for performance under pressure, this customizable platform consistently outperforms in demanding emergency environments—making it the preferred choice for departments that demand vehicles as reliable as their personnel.



Model available as

**SELECT**  
BY BRAUN

# CHIEF XL

## RUGGED & DEPENDABLE

### FULLY WELDED-IN CABINETRY

Crash-tested lightweight cabinets remain secure to protect lives

### SOLID ACRYLIC COUNTERTOPS

High-quality material, enhanced aesthetics, and rounded edges

### ARCHED ROOF DESIGN

Signature roof design adds strength and eliminates water and snow pooling

## FULLY FUNCTIONAL

### EZ GLIDE™ SLIDING DOOR

Enhanced safety and convenience for your crew

### CONFIGURABLE INTERIOR & EXTERIOR CABINETS

Configured cabinets to meet department working style and more

### EBERHARD FREE FLOATING LATCHES

Easier to grip & holds up during everyday use

## BUILT FOR LIFE

### SOLIDBODY™ CONSTRUCTION

Life-saving strength, also efficiently reduces fuel and operational costs

### LIFETIME STRUCTURAL WARRANTY

One of the most comprehensive warranties in the industry, transferable through resale or remount

### REMOUNTABLE

Truly 'Built for Life', ask us about our remount packages



	TYPE I					TYPE III	
	Ford F450*	Ford F550	Chevy 5500HD	RAM R4500	RAM R5500	Ford E450*	Chevy G4500
<b>Fuel Type</b>	Gas or Diesel	Gas or Diesel	Diesel	Gas or Diesel	Gas or Diesel	Gas	Gas
<b>Overall Length</b>	307"	307"	298.75"	306.5"	306.5"	272.75"	279"
<b>Overall Width</b>	98"	98"	98"	98"	98"	98"	98"
<b>Overall Height</b>	112"	112"	112"	112"	112"	112"	112"
<b>Module Length</b>	169"	169"	169"	169"	169"	169"	169"
<b>Module Width**</b>	96"	96"	96"	96"	96"	96"	96"
<b>Headroom</b>	72"	72"	72"	72"	72"	72"	72"
<b>Wheel Base</b>	193"	193"	198"	193"	193"	158"	159"

\*Ford F450 & E450 not available in 4x4 chassis

\*\*Add 2" to overall width to include rub rails & exterior side warning lights



FIND YOUR NEAREST DEALER AT  
**BraunAmbulances.com**



1170 Production Drive  
 Van Wert, OH 45891

Toll Free (877) 344-9990